

ALYS Labware, Lausanne
Tel : 021 312 42 60
Fax : 021 312 42 61
labware@alys-technologies.com



BioPette™ Digital Pipettes

Your research demands the best...

- ❑ *Designed and engineered for maximum accuracy and precision*
- ❑ *Streamlined handle design provides exceptional user comfort*
- ❑ *Adjustable ejector accommodates a wide variety of tips*
- ❑ *Two interchangeable ejector buttons provide personalized comfort*
- ❑ *"Easy-grip" ejector handle for quick removal and adjustment*
- ❑ *Recessed, easy to turn volume adjustment wheel*
- ❑ *7 models cover volume range from 0.1µl to 5000µl*



...Labnet provides it



Interchangeable tip ejector buttons



Ejector removal is a snap with the easy-grip handle



Adjusts to fit various types and brands of tips

Accurate, reliable, comfortable pipetting labware@alys-technologies.com

Your research demands only the highest quality when it comes to pipetting. *BioPette* pipettes meet these quality demands by providing exceptionally accurate and reliable pipetting while standing up to the rigors of everyday use. ISO 9001 certification gives an additional assurance of quality. Their comfortable feel and ergonomic design reduce operator fatigue and down time.

❑ *Precise and accurate*

The pipetting mechanism within each *BioPette* includes a finely polished, stainless steel piston that provides excellent accuracy and reproducibility. The high-precision micrometer is protected from the effects of hand-warming by the pipette's insulated body. To ensure a proper fit for most types and brands of pipette tips, the tip ejector is adjustable. Each pipette is individually calibrated before leaving the factory and will stay in calibration with proper use.

❑ *User friendly*

The streamlined handle of the *BioPette* is comfortable and easy to hold for all users, regardless of hand size. The plunger and tip ejector have been placed to reduce fatigue and repetitive strain injuries. Two interchangeable ejector buttons are supplied to accommodate both large and small hands. The tip ejector is easily removed for pipetting into narrow tubes, yet will not slip during use. At the top of the ejector is a cushioned handle for convenient removal and adjustment of the ejector without straining your fingers. The recessed thumb wheel is easy to turn and recessed to prevent accidental change. Selected volumes are shown on the large, easy to read, digital display.

❑ *Reliable*

BioPettes are solidly constructed to provide years of reliable, trouble-free service. Seals are chemically resistant for a long life. The internal mechanism is self lubricating - no routine maintenance is required. All internal and external parts are constructed of the highest quality materials to withstand the rigors of everyday use.

❑ *Complete range*

Seven models cover the volume range from 0.1µl to 5000µl (the 5000µl model does not have a tip ejector). Each *BioPette* is supplied with a certificate of quality and covered by a one year warranty. For best performance, we recommend the use of *BioPette* pipette tips.



Linear and Carousel
PipetteStands

BioPette Adjustable Pipettes

Cat. No.	Model	Volume Range	Accuracy	Coef. Variation	Tip
P3960-2	BP2	0.1-2µl	±12 to ±1.5%	±6 to ±0.7%	P3984
P3960-10	BP10	0.5-10µl	±3.0 to ±0.5%	±3.0 to ±0.3%	P3984
P3960-20	BP20	0.2-20µl	±4.0 to ±0.8%	±3.0 to ±0.4%	P3990/P3991
P3960-100	BP100	10-100µl	±1.6 to ±0.8%	better than ±1%	P3990/P3991
P3960-200	BP200	20-200µl	±1.5 to ±0.6%	better than ±1%	P3990/P3991
P3960-1000	BP1000	100-1000µl	±0.9 to ±0.6%	better than ±1%	P3995/P3996
P3960-5000*	BP5000	500-5000µl	±0.6 to ±0.5%	better than ±1%	P3998

P3989 Carousel stand for 5 pipettes

P3999 Linear stand for 5 pipettes

*BP5000 does not have a tip ejector

BioPette Tips

P-3984	Pipette tips, natural, 10µl max., 96 tips per rack, 10 racks per pack
P-3990	Pipette tips, yellow, 200µl max., bulk pack of 1000 tips
P-3991	As above, 96 tips/rack, 10 racks with lids
P-3995	Pipette tips, blue, 1000µl max., bulk pack of 1000 tips
P-3996	As above, 100 tips per rack, 10 racks per pack
P-3998	Pipette tips, natural, 5000µl max., bulk pack of 300 tips